SPATZ³

ADJUSTABLE WEIGHT LOSS SYSTEM

EN Instructions for use

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DA Brugsanvisning

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SW Bruksanvisning

Introduction

CAUTION: This device is restricted to sale by or on the order of a licensed physician.

Spatz3 Adjustable Balloon System® Insertion Kit, Cat. No. A-SP3-03K

The **Spatz3 Adjustable Balloon System®** is designed to assist weight loss by partially filling the stomach and inducing satiety.

The **Spatz3 Adjustable Balloon System®** is placed in the stomach and filled with saline, causing it to expand into a spherical shape. The filled balloon is designed to act as a space filler and move freely within the stomach. The unique inflation tube design of the **Spatz3 Adjustable Balloon System®** permits balloon volume adjustment at the time of placement, and at a later time. Balloon volume may range from 400-700 ml (See below for details.)

The **Spatz3 Adjustable Balloon System®** consists of a balloon connected to a curved catheter containing an extendable inflation tube that is connected to an extension tube. The extendable inflation tube exits the catheter behind the balloon, which enables balloon adjustment after the initial implantation.

Indications for Use

The **Spatz3 Adjustable Balloon System®** is indicated for temporary use in weight-loss therapy for overweight and obese patients who have significant health risks related to their weight and have failed to achieve and maintain weight loss with a supervised weight-control program (see specific indications below).

The maximum placement period for the **Spatz3 Adjustable Balloon System®** is 12 months, and must be removed at that time or earlier.

CAUTION: Inflation volumes above 700 ml are not recommended.

Each physician and patient should evaluate the risks associated with endoscopy and intragastric balloons (see **Warnings and Precautions and Complications below**) and the possible benefits of a temporary treatment for weight loss prior to use of the **Spatz3 Adjustable Balloon System®**.

The Spatz3 Adjustable Balloon System® is indicated for:

- Pre-surgical temporary use for weight loss in morbidly obese patients (BMI 40 and above or a BMI of 35 with co-morbidities) prior to bariatric or other surgery, in order to reduce surgical risk.
- Temporary use for weight loss in overweight and obese patients (BMI > 27) who have failed to achieve and maintain weight-loss with a supervised weight control program.

The **Spatz3 Adjustable Balloon System®** is to be used in conjunction with a long-term supervised diet and behavior modification program designed to increase the possibility of long term weight-loss maintenance.

CAUTIONS:

- **Temporary weight-loss treatments have been shown to have poor long-term success rates (weight-loss maintenance) in obese and morbidly obese patients.**

- Physicians have reported the concurrent use of PPI (proton pump inhibitor) and H2 Blocker medications which reduce acid formation or reduce acidity. Silicone elastomer is degraded by acid. Moderating the pH in the stomach with these medications should prolong the integrity of the **Spatz3 Adjustable Balloon System®**, and their use is highly recommended

- **Deflated devices should be removed promptly.** Bowel obstructions have been reported in currently marketed balloons due to deflated balloons passing into the intestines and have required surgical removal. Some obstructions have reportedly been associated with patients who have diabetes or who have had prior abdominal surgery, so this should be considered in assessing the risk of the procedure. Bowel obstructions can result in death.

- The risk of intestinal obstruction may be higher in patients who have had prior abdominal or gynecological surgery. The risk of intestinal obstruction may be higher in patients who have a dysmotility disorder or diabetes.

- Pregnancy or breast-feeding contraindicates use of this device. Should pregnancy be confirmed at any time during the course of treatment, it is recommended that the device be removed.

- Patients should ensure that effective birth control is in place prior to implantation and throughout the duration of implantation, and that they tell their doctor as soon as possible if pregnancy is confirmed during treatment.

- Proper positioning of the **Spatz3 Adjustable Balloon System®** within the stomach is necessary to allow proper inflation. Lodging of the balloon in the esophageal opening during inflation may cause injury and/or esophageal rupture. The physiological response of the patient to the presence of the **Spatz3 Adjustable Balloon System®** may vary depending upon the patient's general condition and the level and type of activity.

- Each patient must be monitored closely during the entire term of treatment in order to detect the development of possible complications. Each patient should be instructed to report to physicians immediately regarding any and all change of symptoms. Symptoms of deflation, gastrointestinal obstruction, ulceration and other complications which might occur should be reviewed with patient, and patients should be advised to contact his/her physician immediately upon the onset of such symptoms.

- Patients with an intragastric balloon that present with severe abdominal pain that have a negative endoscopy and x-ray may additionally require a CT scan to definitively rule out a perforation.

- Patients reporting loss of satiety, increased hunger and/or weight gain should be examined radiographically and/or endoscopically, as this is indicative of a balloon deflation. Should the loss of satiety, increased hunger and/or weight gain not be associated with balloon deflation or other physiological cause, and the balloon is below the recommended maximum (700 ml) capacity, the physician may consider increasing the balloon volume. (See Balloon Volume Adjustment instructions below.)

- If it is necessary to replace a balloon which has spontaneously deflated, the recommended initial fill volume of the replacement balloon is the same as for the first balloon or the most recent volume of the removed balloon. A greater initial fill volume in the replacement balloon may result in severe nausea, vomiting or ulcer formation.

- The **Spatz3 Adjustable Balloon System**[®] is composed of soft silicone elastomer and is easily damaged by instruments or sharp objects. The balloon must be handled only with gloved hands and as recommended in the Instructions for Use.

- The **Spatz3 Adjustable Balloon System**[®] has been tested safely in a 3T MRI scanner. Patients with implanted Spatz3 adjustable balloon systems may safely undergo MRI scanning with a 3T scanner.

Contraindications

- Less than 18 years.
- Treatment represents an unreasonable risk to the patient.
- Pregnancy.
- History of inflammatory disease of the gastrointestinal tract (including ulceration, **any erosive gastritis**, grade 3-4 esophagitis, Barrett's esophagus, or specific inflammation such as Crohn's disease or ulcerative colitis).
- Severe coagulopathy; upper gastro-intestinal bleeding conditions such as esophageal or gastric varices, congenital or acquired intestinal telangiectasia.
- Congenital or acquired anomalies of the GI tract such as atresias or stenoses.
- Emotional disorders, eating disorders, or substance abuse.
- History of a known diagnosis or pre-existing symptoms of systemic lupus erythematosus, scleroderma, or other autoimmune connective tissue disorder.
- Medical illnesses such as, but not limited to cardiopulmonary disease that would increase the risk of endoscopy or weight loss.
- Any other medical condition which could increase the risk of elective endoscopy or weight loss.
- Previous gastric, esophago-gastric, bariatric or duodenal surgery.
- **CLINICALLY SIGNIFICANT HIATAL HERNIA, ESOPHAGEAL DIVERTICULUM, STENOSES, RINGS OR WEBS THAT CAN CAUSE DEVICE IMPLANTATION ABOVE THE DIAPHRAGM.**
- Motility disorders of the GI tract such as gross esophageal motility disorders, gastroparesis or intractable constipation.
- Subjects who require continuous therapy with known ulcerogenic medication (e.g., aspirin, non-steroidal anti-inflammatory agents, COX-2 inhibitors), anti-coagulants or anti-platelet agents.
- Subjects unable or unwilling to participate in a medically supervised diet and behavior modification program with frequent follow up meetings.
- Allergic to Silicone.

Complications

Possible complications of the use of the **Spatz3 Adjustable Balloon System**[®] include, but are not limited to:

- Intestinal obstruction by the balloon. This may require abdominal surgery and/or subsequent surgery. Death due to complications related to intestinal obstruction or surgery is possible.
- Balloon deflation and subsequent replacement.
- Renal Injury

- Esophageal obstruction. Once the balloon has been inflated in the stomach, the balloon could be pushed back into the esophagus. If this occurs, surgery or endoscopic removal could be required. Injury to the digestive tract during placement of the balloon in an improper location such as in the esophagus or duodenum. This could cause bleeding or even perforation, which could require a surgical correction for control.

- Insufficient or no weight loss.
- Injury to teeth, tissue in the oral cavity or throat and upper esophageal sphincter.
- Temporary weight-loss treatments have been shown to have poor long-term success rates (weight-loss maintenance) in severely obese patients.
- Adverse health consequences resulting from weight loss.
- Gastric discomfort, feelings of nausea and vomiting following balloon placement as the digestive system adjusts to the presence of the balloon.
- Continuing nausea and vomiting with dehydration, electrolyte abnormalities, weakness, fainting or falling episode. This could result from direct irritation of the lining of the stomach or as a result of the balloon blocking the outlet of the stomach. It is even theoretically possible that the balloon could prevent vomiting (not nausea or retching) by blocking the inlet to the stomach from the esophagus.
- A feeling of heaviness in the abdomen. Abdominal or back pain, either steady or cyclic.
- Gastroesophageal reflux.
- Pancreatitis (Including Acute).
- Influence on digestion of food. Blockage of food entering into the stomach.
- Bacterial growth in the fluid which fills the balloon. Rapid release of this fluid into the intestine could cause infection, fever, cramps and diarrhea.
- Injury to the lining of the digestive tract as a result of direct contact with the balloon, catheter, polypectomy snare, or as a result of increased acid production by the stomach - esophagitis, gastritis or duodenitis. This could lead to ulcer formation with pain, bleeding or even perforation. Surgery could be necessary to correct this condition, and could result in death.
- Catheter entrance into the duodenum causing anorexia, nausea, vomiting or abdominal pain. This could spontaneously resolve or require endoscopy to pull the catheter back into the stomach, and could become a recurring problem, which would necessitate device removal.
- Spontaneous Hyperinflation is the enlargement of the balloon with extra air that can occur spontaneously. This can lead to symptoms such as pain, nausea, vomiting, dehydration, ulceration, perforation, and could require a down adjustment or removal of the balloon.
- Foul smelling belches or burps.

Complications of routine endoscopy include:

Adverse reaction to sedation or local anaesthetic; Abdominal cramps and discomfort from the air used to distend the stomach; Sore or irritated throat following the procedure; Aspiration of stomach contents into the lungs; Cardiac or respiratory arrest (these are extremely rare and are usually related to

severe underlying medical problems); Digestive tract injury or perforation; Upper digestive tract bleeding.

How Supplied

Each **Spatz3 Adjustable Balloon System**® Insertion Kit contains a balloon attached to an extension tube with a 3-way stopcock at its proximal end and a 50 ml or 60 ml Balloon Syringe. The white cap for closure of the inflation valve is provided within the packaging in a separate pouch. A device insertion facilitator is provided within the packaging.

A patient card that identifies the Spatz3 balloon, is supplied in each package and is to be carried by the patient at all times during the period of balloon implantation.

Spatz3 Adjustable Balloon System adjustment/extraction Kit contains extension tubes with 3-way stopcock and are separately available.

The Spatz3 adjustable balloon Insertion Kit contains: **Spatz3 Adjustable Balloon System**® in a non-sterile state, and the Balloon Syringe in a sterile state

All are supplied **FOR SINGLE USE ONLY** and none should be re-used or sterilized. The risk of infection, incorrect function and patient harm are possible should these devices be re-used or sterilized.

STORE AT ROOM TEMPERATURE (15°-30° C).

Safe Disposal

Dispose of used device/parts in accordance with institution's standard protocol for disposal of biohazard materials

Cleaning Instructions

In the event that the product becomes contaminated prior to use, it should not be used but should be returned to the manufacturer.

CAUTION: DO NOT SOAK THE PRODUCT IN A DISINFECTANT.

The silicone elastomer may absorb some of the solution which could subsequently leach out and cause tissue reaction.

Instructions for use



Use of this device requires a working knowledge of endoscopic techniques and equipment and oral- gastric intubation. The instructions for use of the **Spatz3 Adjustable Balloon System**® do not obviate the need for formal training in these methods. Formal training in the Spatz3 Balloon insertion, down/up adjustment and removal procedures are required and should be performed by a Spatz-authorized instructor. for tutorials movies see Spatzmedical.com

The inflation tube is delicate - please inflate cautiously, since continued inflation while the tube is kinked can result in inflation tube tear. The balloon should not be used if any damage is noted. A stand by **Spatz3 Adjustable Balloon System**® should be available at the time of placement.

Insertion Instructions:

1. Conscious sedation as per standard endoscopy.
2. Perform endoscopy and assess for findings that contraindicate device implantation.
3. Remove the endoscope.

4. Remove the **Spatz3 Adjustable Balloon System**® from the package.
5. Pull the "band-off" to remove **ALL** of the bands.

6. Dry the endoscope thoroughly.
7. Place the insertion facilitator on the distal end of the endoscope and roll it over the distal edge of the balloon. The tip of the balloon should be within 3 mm of the proximal hole of the insertion facilitator-as shown in photos. Push the Valve Hold snugly into the cap of the white catheter (see figures below)

WARNING-Skipping this step can result in an inflation tube tear which will require removal of the balloon. WARNING: Do not twist the valve during this step it may cause a kink in the inflation tube and failure to inflate.

Connect the valve firmly to the extension tube. Lubricate the device generously over the length of the balloon, and the inflation valve, using medical grade jelly or equivalent. Pass the endoscope and device in one hand as one unit, into the throat and advance to the esophagus and stomach as per usual endoscopic insertion. Once in the gastric antrum, retroflex and confirm that the balloon has passed into the stomach, connect one opening of the 3-way stopcock to the large syringe and a second opening to the solution bag and inflate the balloon with 500 ml. During inflation, if there is any resistance to inflation, reposition the endoscope until the inflation flows easily.

1. place condom

2. pull proximally


3. pull the other side


4. until edge is at the tip of scope


5. Roll back to cover proximal hole (edge may go beyond scope tip)


6. Balloon tip touches condom - begin rolling

7. Roll back fully (balloon tip must be within 3 mm of proximal hole)

8. Pull proximal hole over scope and then distal hole

9. Push ends back toward balloon tip

10. Valve-Hold

11. Push Valve-Hold into white cup (push straight – do not twist)

12. Valve-Hold inside cup

13. Ready for use

8. Inflate or deflate to the final desired volume—we recommend starting with 450-500 ml. Then remove the endoscope fully – you may feel some mild pressure as the balloon separates from the Insertion Facilitator.
9. Pull out the extension tube until the valve exits the mouth. The valve is grasped firmly with a gauze pad while the extension tube is twisted off and replaced with the white cap. The white cap is provided in a separate bag in the packaging.
Warning: The cap and valve are slippery, and it is difficult to be sure that you have closed the cap fully. You must dry the valve and cap with paper towel before closing the cap.
10. Release the capped valve into the back of the throat.
11. The white cap and white catheter is pushed down below the GE junction by the tip of the endoscope. Remove the endoscope.

Balloon Volume Adjustment instructions:

1. An extension tube with 3 way stopcock is provided in the “Spatz Adjustable Balloon System® Adjustment/ Extraction Kit,” Cat. No. A-SP3-015, or can be self supplied.
2. Perform endoscopy and inspect for balloon or inflation tube leakage.
3. Capture the valve loop with a rat tooth forceps.
4. Remove the endoscope with valve pulled against the endoscope tip. **THE VALVE MUST BE TOUCHING THE ENDOSCOPE TIP DURING REMOVAL TO PREVENT IT FROM GETTING CAUGHT AT THE GE-JUNCTION.** Once the tip of the endoscope with the ensnared valve exits the mouth, grab the valve firmly. One should feel the elastic recoil pressure of the stretched inflation tube. This step may require another set of hands.
5. Twist off the white cap from the valve and attach the valve to the provided male port of the extension tube and securely close it **DO NOT LOSE THE WHITE CAP!**
6. The extension tube is slowly allowed to slide into the throat, until the inflation tube is no longer taut.
7. **REINSERT THE ENDOSCOPE TO CONFIRM POSITION OF THE BALLOON IN THE STOMACH.**
8. Attach the 50-60 ml Syringe (self supplied) and the solution bag to the 3-way stopcock and commence inflation/deflation (range 400-700 ml).
9. **NEVER INFLATE/DEFLATE WITHOUT ENDOSCOPIC VIEW OF THE INFLATION TUBE.** If the inflation tube is kinked or folded, use the tip of the endoscope or movement of the extension tube to unkink the inflation tube.
10. When adjustment is completed, leave the endoscope in the mid-esophagus while removing the extension tube fully until the valve exits the mouth. The valve is grasped firmly while the extension tube is twisted off and replaced with the white cap.
Warning: The cap and valve are slippery, and it is difficult to be sure that you have closed the cap fully. You must dry the valve and cap with paper towel before closing the cap.

- The capped valve is lowered into the back of the throat and allowed to glide down the esophagus.
- Confirm the passage of the valve into the esophagus. It can get caught in the mouth between the teeth and the tongue -if so a finger sweep should release it. The white cap and white catheter is pushed down below the GE junction. Remove the endoscope.

Device Removal Instructions:

- An extension tube with 3 way stopcock is provided in the "Spatz Adjustable Balloon System® Adjustment/Extraction Kit," Cat. No. A-SP3-015, or can be self supplied. Or a needle catheter for balloon puncture can be used.
- Perform endoscopy, at which time the balloon can be deflated by needle puncture or utilizing the inflation tube (steps 3 to 8 in the adjustment section above, and disconnect extension tube before balloon extraction) or both.
- DEFLATION USING THE INFLATION TUBE REQUIRES ENDOSCOPIC OBSERVATION. (avoid tube kinking)**
- Allow the balloon to completely empty of all the fluid.
- A grasping forceps or snare may be used to extract the balloon.
- Pull the balloon up against the tip of the endoscope, and remove the endoscope. Hyperextension of the neck of the patient can help with removal as the device passes through the upper esophageal sphincter.

Warranty and Limitation of Warranty

Spatz FGIA, Inc. warrants that reasonable care was used in the design and manufacture of this product. Because Spatz FGIA, Inc. has no control over the conditions of use, patient selection or handling of the device after it leaves its possession, Spatz FGIA, Inc. does not warrant either a good effect or against an ill effect following its use. Spatz FGIA, Inc. shall not be directly or indirectly responsible for any incidental or consequential loss, damage or expenses directly or indirectly arising from the use of this product. Spatz FGIA, Inc.'s SOLE responsibility in the event Spatz FGIA, Inc. determines the product was defective when shipped by Spatz FGIA, Inc., shall be the replacement of the product. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including but not limited to any implied warranties of merchantability or fitness for use.

Product Specifications

Spatz3 Adjustable Balloon System® Insertion Kit, Catalog No. A-SP3-03K:

Balloon, cap, extension tube with 3-way stopcock insertion facilitator (Note: The products are supplied non-sterile and packaged for use), and Balloon syringe (sterile)

Spatz3 Adjustable Balloon System® Adjustment/Extraction Kit, Cat. No. A-SP3-015 contains extension tubes with 3-way stopcock.

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For additional information, please contact your local distributor or:


Spatz FGIA Inc.
1801 S Perimeter Rd, Suite 130
Fort Lauderdale, FL 33309 USA

Tel: 1-516-303-0613
Fax: 1-516-303-0614
www.spatzmedical.com

MedNet EC-REP IIb GmbH
Borkstrasse 10, 48163
Münster, Germany

Label Symbol Glossary

REF	Device reference number	EC REP	Authorized Representative in the European Community	LATEX FREE	Device does not contain natural latex rubber
LOT	Lot number in the format of LLLLLL	!	Attention, see instructions for use	Rx only	The device is restricted to sale by or on the order of a physician
Manufacturing date icon	Manufacturing date in the format of YYYY-MM	NON STERILE	Non-sterile	STERILE EO	sterilized using ethylene oxide
Expiry date icon	2-year expiry date, based on the manufacturing date, in the format of YYYY-MM	CE 1023	CE marking number of the notified body	Limits of temperature icon	Limits of temperature
Manufacturer icon	Manufacturer	Do not reuse icon	Do not reuse	Consult instruction of use icon	Consult instruction of use
MD	Medical Device	UDI	Unique Device Identifier	Don't use if package is damaged icon	Don't use if package is damaged
Health care center or doctor icon	Health care center or doctor	Implantation Date icon	Implantation Date	Patient identification icon	Patient identification
				Patient information website icon	Patient information website